




Schutz nach DIN 24 beansprucht.
Use of this document without permission
without express authorisation is forbidden.
Dräger Medical Technik GmbH, Lübeck



- Nur aufkleben in Verbindung mit Hubsaeule 2M20946
- ⊗ Nur aufkleben bei BT 8000 WB mit Matr.-Heizung oder bei BT 8004/8010
- ◐ Nur aufkleben bei BT 8000 WB
- ◑ Nur aufkleben bei BT 8000 mit fester Saeule
- ⑦ ⊗ ⊗ Nur aufkleben, wenn nicht 2M22263 aufgeklebt ist.

* enthalten in
Schildersatz Babytherm 8000 WB USA
Babytherm 8004/8010 USA 2M21896 ST

scale	general tolerances	material	safety class
		SM-Papier DIN A2 ws hf 80g/qm	<e> 0 <m> 0
field reserved		revision text	
		Darstellung geändert, Zeichn.-Nr., Masz und Text neu; Zeichn.-Nr. entfallen, Masz geändert.	(2000 / CC / 289)
filing dept. 2950	responsible dept. 5511	author 23.11.2000 Stickel	approved 23.11.2000 Hartwig
type drawing	language de	title Klebeanweisung	contract n
	replacement for 2M21901-06	number 2M21901	
creating system CADRA	release status RELEASED	revision index 07	occurs 4x
		page/of 01/01	Dräger Medizintechnik GmbH D-23542 Lübeck
			size A2

CAUTION!

KEEP SPACE
BELOW
CABINET
CLEAR WHEN
LOWERING

KEEP SPACE
ABOVE
BASSINET
CLEAR WHEN
RAISING

CAUTION!

KEEP SPACE
BELOW
CABINET
CLEAR WHEN
LOWERING

KEEP SPACE
ABOVE
BASSINET
CLEAR WHEN
RAISING

DANGER!

RISK OF EXPLOSION IF USED IN THE
PRESENCE OF FLAMMABLE ANESTHETICS

WARNING!

DISCONNECT SUPPLY BEFORE SERVICING

REPAIRS ON THIS EQUIPMENT TO BE PERFORMED
ONLY BY DrägerService OR ITS AUTHORIZED
SERVICE CENTERS

CAUTION!

TO MAINTAIN GROUNDING INTEGRITY, CONNECT
ONLY TO A "HOSPITAL GRADE" RECEPTACLE
TO REDUCE RISK OF ELECTRIC SHOCK, DO
NOT REMOVE COVER

FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO
SALE BY OR ON THE ORDER OF A PHYSICIAN

FOR ACTIVE AND BIGGER BABIES
USE SIDE PANELS OF 23 CM HEIGHT
AND ADDITIONAL INNER PANELS.
MAKE SURE TO MONITOR MORE ACTIVE
BABIES APPROPRIATELY

WARNING!

DISCONNECT BOTH SUPPLIES
BEFORE SERVICING

WARNING!

DISCONNECT BOTH SUPPLIES
BEFORE SERVICING

WARNING!

DANGER OF PATIENT INJURY
NEVER LEAVE BABY UNATTENDED WHEN
SURROUNDING PANELS ARE DOWN

WARNING!

DANGER OF PATIENT INJURY
NEVER LEAVE BABY UNATTENDED WHEN
SURROUNDING PANELS ARE DOWN

WARNING!

WHEN USED WITH THE RADIANT WARMER CHECK ITS
OPERATORS MANUAL FOR DIRECTIONS OF THE
CORRECT HEAT LEVEL AND DISTANCE TO THE BABY

WARNING!

DO NOT SWITCH OFF BABY THERM WHILE A BABY IS LYING ON THE
GEL MATTRESS AS THIS WILL RESULT IN COOLING OF THE PATIENT

WARNING!

CHECK BABY'S TEMPERATURE AND SKIN
CONDITION AT REGULAR INTERVALS AND ADJUST
TEMPERATURE SETTING TO INDIVIDUAL PATIENT NEEDS

WARNING!

USE FOAM RUBBER MATTRESS WHEN THE HEATING
MATTRESS IS COOL - GEL MATTRESS WILL COOL THE INFANT IF
MATTRESS HEATER IS SWITCHED OFF

WARNING!

FOR ACTIVE AND BIGGER BABIES USE SIDE PANELS
OF 23 CM HEIGHT AND ADDITIONAL INNER PANELS
MAKE SURE TO MONITOR ACTIVE BABIES APPROPRIATELY

DRÄGER MEDICAL, INC.
3135 QUARRY ROAD
TELFORD, PA 18969
TEL: (215) 721-6910
FAX: (215) 721-6915

WARNING!

DO NOT TOUCH OR COVER TOP OF RADIANT HEATER
DO NOT TOUCH PROTECTIVE SCREEN.
DANGER OF BURN INJURIES!

WARNING!


DO NOT TOUCH OR COVER TOP OF RADIANT HEATER
DO NOT TOUCH PROTECTIVE SCREEN.
DANGER OF BURN INJURIES!

Schildersatz Babytherm 8000 WB USA 2M21896
Babytherm 8004/8010 USA 2M21896

- auf einem Bogen geliefert
Abmessung und Anordnung sind dem
Hersteller ueberlassen
- Folie zwischen den Schildern
entfernt

Hersteller Fa. Seridruck

ST2M21896



2M21898(2x)

2M21894

scale	1:1	general tolerances	material	safety class
field reserved			revision text	<e> 0 <m> 0
			Text geändert	
filing dept.	2950	responsible dept.	5511	(2000 / CC / 409)
author	01.12.2000 Stickel	approved	01.12.2000 Hartwig	contract
date				n
name				
title				
type	drawing	language	de	Schildersatz
replacement for	2M21896-05	number	2M21896	
creating system	CADRA	release status	RELEASED	
revision index	06	occurs	1x	page/of
				01/01
				size A2

Dräger

Dräger Medizintechnik
GmbH
D-23542 Lübeck